

## Cytoprep-1 Tabletop Cyto Liquid Smear Centrifuge

Ideal for TCT thin layer cytology test and LBC liquid based cytology test.

### Features :

- Clear blue tempered glass top lid, safe & compact, it consists of microprocessor control
- Brush-less DC motor, free maintenance, no powder pollution
- LCD display Speed/RCF & Time
- There are 10 rates of speed up and down faster and stable.
- Safety door lock alarm inside
- There are five angle rotors are for options
- ISO9001:2008, ISO13485, & CE certificates are available as well as calibration report.




### Technical Data

Max.Speed	4,000r/min
Max.Capacity	6x9ml
Timer	0~99min
Dimension	540x370x280mm
Speed Accuracy	±20r/min
Max.RCF	2,170xg
Power supply	AC220V, 50/60Hz
Noise	≤55dBA
Net /Gross weight	19.5/21.5KG
Temperature range	Non-refrigerated



Adapter for 15ml conical tube



### Matched Rotors

Order No.	Rotor type	Max.Speed (rpm)	Max.Volume (ml)	Max.RCF (xg)
No21427	Cyto rotor	4,000	6x9ml	2,170xg
	Adapter matching swing out rotor body	4,000	6x15ml	2,250xg

