

TD5 Tabletop Low Speed Centrifuge

The capacity is 4 x 250ml in Swing out rotor matching Adapters for 5ml, 7ml, 10ml, 50ml tubes...to meet various of testings flexible. Widely used in the fields of biochemistry, biological products, pharmaceutical factory and laboratory.

Features :

- metal housing is safe and compact. It is of microprocessor control accuracy high.
- Brush-less frequency motor, free maintenance, no powder pollution
- Digital display Speed/RCF & Time for options
- There are 10 rates of speed up and down faster and stable.
- Electronic safety lock alarm inside, just press down the top lid until it will lock in sealed and open up after separation automatically.
- having both swing out rotors and angle rotors are for options
- ISO9001:2008, ISO13485, CE certificates are available.

Technical Data:



Max.Speed	5,000r/min
Max.Capacity	4x250ml
Timer	1~99min
Dimension	570x400x320mm
Speed Accuracy	±20r/min
Max.RCF	4,730xg
Power supply	AC220V, 50/60Hz
Noise	≤58dBA
Net weight	43KG
Temperature range	Non-refrigerated



Matched Rotors

Order No.	Rotor type	Max.Speed (rpm)	Max.Volume (ml)	Max.RCF (xg)
No30569	Swing out rotor	4,000	4x250ml	2,880xg
Adapters in 250ml buckets (round cup)				
For 5ml tube to do 4x7x5ml	For 7ml tube to do 4x7x7ml	For 10ml tube to do 4x7x10ml	For 15ml tube to do 4x4x15ml	For 50ml tube to do 4x1x50ml
No31491	Microlater rotor	4,000	2x4x96well	2,910xg
No31377	Swing out rotor	5,000	4x50ml	4,730xg
No31378		4,000	4x100ml	4,730xg
		4,000	32x10ml	2,580xg
		4,000	24x15ml	3,020xg
			64x5ml	2,470xg
			64x7ml	2,860xg
No30627	Angle rotor	5,000	30x15ml	3,830xg
No30641	Angle rotor	5,000	12x50ml	3,860xg
No30614	Angle rotor	5,000	6x100ml	3,130xg