


TXK4 Low Speed Centrifuge For Blood Gel Cards

Widely used in the fields of blood type serology., blood type testing, Micro column gel & immune detection etc.,

Features :

- Metal housing is safe and compact. It is of microprocessor control accuracy high.
- Brush-less DC motor, free maintenance, no powder pollution
- Microprocessor control and LCD display Speed/RCF & Time
- There are 10 rates of speed up and down faster and stable.
- Safety door lock alarm inside,
- There are five angle rotors are for options
- ISO9001:2008, ISO13485, CE certificates are available as well as calibration report.



Technical Data (TXK4-12cards)

Max.Speed	4,000r/min
Max.Capacity	12cards
Timer	0~99min
Dimension	485x320x250mm
Speed Accuracy	±20r/min
Max.RCF	1,480xg
Power supply	AC220V, 50/60Hz
Noise	≤55dBA
Net weight	17KG
Temperature range	Non-refrigerated

TXK4-24 cards

Max.Speed	4,000r/min
Max.Capacity	24cards
Timer	0~99min
Dimension	500x360x250mm
Speed Accuracy	±20r/min
Max.RCF	2,120xg
Power supply	AC220V, 50/60Hz
Noise	≤55dBA
Net weight	25KG
Temperature range	Non-refrigerated



Matched Rotors

Order No.	Rotor type	Max.Speed (rpm)	Max.Volume(ml)	Max.RCF (xg)
No33183	Angle rotor	3,000	12cards	830xg
No33184	Angle rotor	3,000	24cards	1,080xg