

## iSeal-200 vertical Intelligent Tube sealer

It is mainly used in blood bank stations, hospitals, pharmaceutical factories, laboratories, disease control centers, etc.,



### Triple Heat Dissipation

Continuous Heat Sealing Capability



### Wide Sealing Surface

Maximum Sealing Tube Diameter



### No Preheating Required

Ready to use upon start up



### Splash Protection Plate

Metal shielding casing for safety



### Can Be Used in Parallel with Multiple Units

Ultra-Slim body design

## Technical Parameters:

Heat Sealing Time	0.5~2s	Heat Sealing tube (Outer diameter)	2~6mm
Indicator Light Colors	Red/Blue	RF Frequency	40.68MHz
Grounding Protection	Class I	Spacing Adjustment	65,70,75,80,85,90,95,100mm
Power Consumption	250w during operation,10W in standby	Power supply	100V~120V/220V~240V AC
Weight	4.5kg	Dimensions	336×65×173mm
Operating Temperature	0~40℃	Storage temperature	-20~70℃